

Recombinant Human C-X-C Motif Chemokine 1,CXCL1 (C-6His)



Catalog No: #AP75277

Orders: order@signalwayantibody.com

Package Size: #AP75277-1 10ug #AP75277-2 50ug #AP75277-3 500ug #AP75277-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human C-X-C Motif Chemokine 1,CXCL1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09341
Uniprot	P09341
GeneID	2919;
Calculated MW	8.9kDa
Target Sequence	ASVATELRCQCLQTLQGIHPKNIQSVNVKSPGPHCAQTEVIATLKNRKAACLNPAPIVKKIIEKMLNSDKSNVD HHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 5% Trehalose, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only