

Recombinant Human C-C motif chemokine 17,CCL17 (C-6His)

Catalog No: #AP75278

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75278-1 10ug #AP75278-2 50ug #AP75278-3 500ug #AP75278-4 1mg

Description

Product Name	Recombinant Human C-C motif chemokine 17,CCL17 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q92583
Uniprot	Q92583
GeneID	6361;
Calculated MW	9.1kDa
Target Sequence	ARGTNVGRECCLEYFKGAIPLRKLKTWYQTSEDCSRDAIVFVTVQGRAICSDPNNKRVKNAVKYLQSLERSVD HHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only