

Recombinant Human Endothelin-converting Enzyme 1,ECE-1 (N-8His)



Catalog No: #AP75279

Orders: order@signalwayantibody.com

Package Size: #AP75279-1 10ug #AP75279-2 50ug #AP75279-3 500ug #AP75279-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Endothelin-converting Enzyme 1,ECE-1 (N-8His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P42892
Uniprot	P42892
GeneID	1889;
Calculated MW	78.8kDa
Target Sequence	HHHHHHHHQYQTRSPSVCLSEACVSVTSSILSSMDPTVDPCHDFFSYACGGWIKANVPDGHSRWGTFSNL WEHNQAIKHLLLENSTASVSEAERKAQVYRACMNTRIEELRAKPLMELIERLGGWNITGPWAKDNFQDTLQ VVTAHYRTSPFFSVYVSADSKNSNSNVIQVDQSGGLPSRDYILNKTENEKVLTYLNYMVQLGKLLGGGDE EAI RPQMQQILDFETALANITIPQEKRRDEELIYHKVTAELQTLAPAINWLPFLNTIFYPVEINESEPIVVYDKEYL EQISTLINTDRCLLNMYMIWNLVRKTSSFLDQRFQDADEKFMVEMYGTTKCLPRWKFCVSDTENNLGFALG PMFVKATFAEDSKSIATEIILEIKKAFEESSLTKWMDDEETRKSAREKADAIYNMIGYPNFMIDPKELDKVFN DYT AVPDL YFENAMRFFNFSWRVTADQLRKAPNRDQWSMTPPMVNAYYSPTKNEIVFPAGILQAPFYTRSSPKAL NFGGIGVVVGHETHAFDDQGREYDKDGNLRPWWKNSSVEAFKRQTECMVEQYSNYSVNGEPVNGRHTLG ENIADNGGLKAA YRAYQNWWKNGAEHSLPTLGLTNNQLFFLGFAQVWCSVRTPESSHEGLITDPHSPSRFR VIGLSNSKEFSEHFRCPPGSPMNPPHKCEVW
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only