

Recombinant Human EDIL3,Del-1 (C-6His)

Catalog No: #AP75280



Package Size: #AP75280-1 10ug #AP75280-2 50ug #AP75280-3 500ug #AP75280-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EDIL3,Del-1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8N610
Uniprot	Q8N610
Calculated MW	53.09kDa
Target Sequence	VPQFGKGDICDPNCPENGGICLPGLADGFSFCECPDGFTDPNCSSVVEVASDEEEPTSAGPCTPNPCHNGG TCEISEAYRGDTFIGYVCKCPRGFNGIHCQHNINECEVEPCKNGGICTDLVANYSCPCGEGFMGRNCQYKCS GPLGIEGGIISNQQTASSTHRALFGLQKWYPYARLNKKGLINAWTAAENDRWPWIQINLQRKMRVTGVITQG AKRIGSPEYIKSYKIAYSNDGKTWAMYKVKGTNEDMVFRGNIDNNTPYANSFTPIKAQYVRLYPQVCRRHCT LRMELLGCELSGCSEPLGMKSGHIQDYQITASSIFRTLNMDFMFTWEPRKARLDKQKGVNAWTSGHNDQSQW LQVDLLVPTKVTGIIQTGAKDFGHVQFVGSYKLAYSNDGEHWTVYQDEKQRKDKVFQGNFDNDTHRKNVIDP PIYARHIRLPWSWYGRITLRSELLGCTEEEDHDDHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only