

Recombinant Human Fc Receptor-Like Protein 1, FCRL1, CD307a (C-6His)

Catalog No: #AP75281

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75281-1 10ug #AP75281-2 50ug #AP75281-3 500ug #AP75281-4 1mg

Description

Product Name	Recombinant Human Fc Receptor-Like Protein 1, FCRL1, CD307a (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96LA6
Uniprot	Q96LA6
GeneID	115350;
Calculated MW	32.24kDa
Target Sequence	AELFLIASPSHPTEGSPVTLTCKMPFLQSSDAQFQCFFRDTRALGPGWSSSPKLQIAAMWKEDTGSYWCEA QTMASKVLRSRRSQINVHRVPVADVLETQPPGGQVMEGDRLVLICSVAMGTGDITFLWYKGAVGLNLQSKT QRSLTAEYEIPSVRESDAEQYYCVAENGYGSPSGLVSITVRIPVSRPILMLRAPRAQAAVEDVLELHCEALRG SPPILYWFYHEDITLGSRSAPSGGGASFNLSLTHEHSGNYSCANNGLGAQRSEAVTLNFTVPTGARSNVDH HHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only