

Recombinant Human Glioma Pathogenesis-Related Protein 1, GLIPR1, RTVP1 (C-6His)



Catalog No: #AP75282

Orders: order@signalwayantibody.com

Package Size: #AP75282-1 10ug #AP75282-2 50ug #AP75282-3 500ug #AP75282-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Glioma Pathogenesis-Related Protein 1, GLIPR1, RTVP1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P48060
Uniprot	P48060
GeneID	11010;
Calculated MW	25.16kDa
Target Sequence	ANILPDIENEDFIKDCVRIHNKFRSEVKPTASDMLYMTWDPALAQIAKAWASNCQFSHNTRLKPPHKLHPNFTS LGENIWTGSPVIFSVSSAITNWyDEIQDYDFKTRICKKVCGHYQVWWADSYKVGCAVQFCPKVSGFDALSNG AHFICNYGPGGNYPTWPYKRGATCSACPNNDKCLDNLCVNRQRDQVKRYYSVVYPGWPIYPRNRVDHHHH HH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only