

Recombinant Human IL-20 Receptor Subunit α,IL20RA (C-6His)

Catalog No: #AP75283

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75283-1 10ug #AP75283-2 50ug #AP75283-3 500ug #AP75283-4 1mg

Description

Product Name	Recombinant Human IL-20 Receptor Subunit α,IL20RA (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UHF4
Uniprot	Q9UHF4
GeneID	53832;
Calculated MW	26.3kDa
Target Sequence	VPCVSGGLPKPANITFLSINMKNVLQWTPPEGLQGKVKVITYTVQYFIYGQKKWLNKSECRNINRTYCDLSAETS DYEHQYYAKVKAIWGTKCSKWAESGRFYPFLETQIGPPEVALTTDEKSISVVLTAPEKWKRNPEDLPVSMQQI YSNLKYNVSVLNTKSNRTWSQCVTNHTLVLTWLEPNTLYCVHVESFVPGPPRAQPSEKQCARTLKDQSSEF KAKVDHHHHHH
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only