

# Recombinant Human Leucine-Rich Repeat-Containing Protein 2, LRRN2 (C-6His)



Catalog No: #AP75284

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75284-1 10ug #AP75284-2 50ug #AP75284-3 500ug #AP75284-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Leucine-Rich Repeat-Containing Protein 2, LRRN2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O43300
Uniprot	O43300
GeneID	26045;
Calculated MW	45.55kDa
Target Sequence	CPPKCRCEKLLFYCDSQGFHSPNATDKGSLGLSLRHNHITELERDQFASFSQLTWLHLDHNQISTVKEDAFQ GLYKLELILSSNKIFYLPNTTFTQLINLQNLDSLNFQSSSLHPELFYGLRKLQTLHLRSNSLRITIPVRLFWDCRS LEFLDLSTNRLRSLARNGFAGLIKREHLEHNQLTKINFAHFLRLSSLHTLFLQWNKISNLTCGMEWTWGTL KDLDTGNEIKAIDLTVFETMPNLKILLMDNKNLNSLDSKILNSLRSLLTVGLSGNLWECSARICALASWLGSFQG RWEHSILCHSPDHTQGEDILDAVHGFQLCWNLSTTVTMATTYRDPTEYTKRISSSSYHVGDKEIPTTAGIAV TTEEHFPEPDNAIFTQRVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only