

Recombinant Human LYVE-1,HAR,XLKD1 (C-6His)

Catalog No: #AP75285



Package Size: #AP75285-1 10ug #AP75285-2 50ug #AP75285-3 500ug #AP75285-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LYVE-1,HAR,XLKD1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y5Y7
Uniprot	Q9Y5Y7
GeneID	10894;
Calculated MW	24.6kDa
Target Sequence	LVQGSRAEELSIQVSCRIMGITLVSKKANQQLNFTEAKEACRLLGLSLAGKDQVETALKASFETCSYGWVGD GFVVISRISPNPKCGKNGVGLIRKVPVSRQFAAYCYNSSDTWTNSCIPEIITTKDPIFNTQTATQTTEFIVSDST YSVASPYSTIPAPTTTTPPAPASTSIPRRKKLICVTEVFMETSTMSTETEPFVENKAAFKNEAAGFGGVPTVDHH HHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris-Citrate, 1 50mM NaCl, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only