

Recombinant Human PAF-AH,PLA2G7 (C-6His)

Catalog No: #AP75286



Package Size: #AP75286-1 10ug #AP75286-2 50ug #AP75286-3 500ug #AP75286-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PAF-AH,PLA2G7 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q13093
Uniprot	Q13093
GeneID	7941;
Calculated MW	48.8kDa
Target Sequence	FDWQYINPVAHMKSSAWVNKIQVLMAAASFGQTKIPRGNPYSVGCTDLMFDHTNKGTFRLRYPSQDNDRL DTLWIPNKEYFWGLSKFLGTHWLMGNILRLLFGSMTPANWNSPLRPGEKYPLVVFHSHGLGAFRTLYSAIGIDL ASHGFIVAAVEHRDRSASATYFQKQSAEIGDKSWLYLRTLKQEEETHIRNEQVRQRAKECSQALSLILDIDH GKPVKNALDLKFDMEQLKDSIDREKIAVIGHSFGGATVIQTLSEDQRFRCGIALDAWMFPLGDEVYSRIPQPLF FINSEYFQYPANIIKMKKCYSPDKERKMITIRGSVHQNFADFTFATGKIIGHMLKKGDIDSNAIDLSNKASLAF QKHLGLHKDFDQWDCLIEGDENLIPGTNINTTNQHIMLQNSSGIEKYNVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM HAc-NaCl, 150mM NaCl, 10% Glycerol, pH 4.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only