

Recombinant Human Proteinase 3, PRTN3, Myeloblastin (C-6His)

Catalog No: #AP75287

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75287-1 10ug #AP75287-2 50ug #AP75287-3 500ug #AP75287-4 1mg

Description

Product Name	Recombinant Human Proteinase 3, PRTN3, Myeloblastin (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P24158
Uniprot	P24158
GeneID	5657;
Calculated MW	25.64kDa
Target Sequence	AEIVGGHEAQPMSRPYMASLQMRGNPGSHFCGGTLIHPSFVLTAHCLRDIPQRLVNVVLGAHNVRTQEPTQ QHFSVAQVFLNNYDAENKLNVDVLLIQLSSPANLSASVATVQLPQQDQPVPHGTQCLAMGWGRVGAHDPPAQ VLQELNVTVVTFFCRPHNICTFVPRRKAGICFGDSSGPLICDGIQIGIDSFVIWGCATRLFPDFTRVALYVDWIR STLRRVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 10mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only