

Recombinant Human Serum Amyloid A1,SAA1 (N-6His)

Catalog No: #AP75290

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75290-1 10ug #AP75290-2 50ug #AP75290-3 500ug #AP75290-4 1mg

Description

Product Name	Recombinant Human Serum Amyloid A1,SAA1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P0DJ18
Uniprot	P0DJ18
GeneID	6288;
Calculated MW	13.2kDa
Target Sequence	MNHKVVHHHHHMRFFSFLGEAFDGARDMWRAYSDMREANYIGSDKYFHARGNYDAAKRGPGGWAAEAI SDARENIQRFFGHGAEDSLADQAADEWGRSGKDPNHFRPAGLPEKY
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM EDTA, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only