

Recombinant Human Syndecan-2,SDC2,CD362 (C-6His)

Catalog No: #AP75291

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75291-1 10ug #AP75291-2 50ug #AP75291-3 500ug #AP75291-4 1mg

Description

Product Name	Recombinant Human Syndecan-2,SDC2,CD362 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P34741
Uniprot	P34741
GeneID	6383;
Calculated MW	14.98kDa
Target Sequence	ESRAELTSDKMYLDNSSIEEASGVYPIDDDYASASGSGADEDVESPELTTTRPLPKILLTSAAPKVETTTLNI QNKIPAQTKSPEETDKEKVLHLSDSERKMDPAEEDTNVYTEKHSDSLFKRTEVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris-Citrate, 150mM NaCl, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only