

Recombinant Human Platelet Endothelial Cell Adhesion Molecule, PECAM-1, CD31 (C-Fc)



Catalog No: #AP75292

Orders: order@signalwayantibody.com

Package Size: #AP75292-1 10ug #AP75292-2 50ug #AP75292-3 500ug #AP75292-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Platelet Endothelial Cell Adhesion Molecule, PECAM-1, CD31 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16284
Uniprot	P16284
GeneID	5175;
Calculated MW	91.6kDa
Target Sequence	QENSFTINSVDMKSLPDWTVQNGKNLTLQCFADVSTTSHVKPQHMLFYKDDVLFYNISSMKSTESYFIPEVRI YDSGTYKCTVIVNNEKTTAEYQVLVEGVSPRVTLDKKEAIQGGIVRVNCSVPEEKAPIHFTIEKLELNEKMMVK LKREKNSRDQNFVILEFPVEEQDRVLSFRCQARIISGIHMQTSESTKSELVTVTESFSTPKFHISPTGMIMEGAQ LHIKCTIQVTHLAQEFPEIIIIQKDKAIVAHNRHGKAVYSVMAMVEHSGNYTCKVESSRISKVSSIVVNITELFSKP ELESSFTHLDQGERLNLSCSIPGAPPANFTIQKEDTIVSQTQDFTKIASKSDSGTYICTAGIDKVVKKSNTVQIVV CEMLSQPRISYDAQFEVIKQGTIEVRCEISISGLTPISYQLLKTSKVLENSTKNSNDPAVFKDNPTEDVEYQCVAD NCHSHAKMLSEVLRVKVIAPVDEVQISILSSKVVESGEDIVLQCAVNEGSGPITYKFYREKEGKPFYQMTSNAT QAFWTKQKANKEQEGEYYCTAFNRANHASSVPRSKILTVRVILAPWKKVDDIEGRMDEPKSCDKTHTCPPCP APELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRV VSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPS DIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPG K
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only