

Recombinant Human Gastric Triacylglycerol Lipase,LIPF (C-6His)

Catalog No: #AP75300

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75300-1 10ug #AP75300-2 50ug #AP75300-3 500ug #AP75300-4 1mg

Description

Product Name	Recombinant Human Gastric Triacylglycerol Lipase,LIPF (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07098
Uniprot	P07098
GeneID	8513;
Calculated MW	44.22kDa
Target Sequence	LFGKLHPGSPEVTMNISQMITYWGYPNEEYEVVTEDEGYILEVNRIPYGKKNSGNTGQRPVFLQHGLLASATN WISNLPNNSLAFILADAGYDVWLGNSRGNTWARRNLYSPDSVEFWAFSFDMAKYDLPATIDFIVKKAGQKQ LHYVGHSQGTIGFIAFSTNPSLAKRIKTFYALAPVATVKYTKSLINKLRFVPSLFFKIFGDKIFYPHNFFDQFLA TEVCSREMLNLLCSNALFIICGFDKFNFTSRLDVYLSHNPAGTSVQNMFWHTQAVKSGKFQAYDWGSPVQN RMHYDQSQPPYYNVTAMNVPIAVWNGGKDLLADPQDVGLLLPKLPNLIYHKEIPFYNHLDFIWAMDAPQEVYN DIVSMISEDKKVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 25mM TrisHCl, 100mM glycine, 10%glycerol,pH 7.3.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only