

Recombinant Human Platelet-derived Growth Factor Receptor Alpha, PDGF R α (C-6His)



Catalog No: #AP75301

Orders: order@signalwayantibody.com

Package Size: #AP75301-1 10ug #AP75301-2 50ug #AP75301-3 500ug #AP75301-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Platelet-derived Growth Factor Receptor Alpha, PDGF R α (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16234
Uniprot	P16234
GeneID	5156;
Calculated MW	57kDa
Target Sequence	QLSLPSILPNENEKVVQLNSSFLRCFGESEVSWQYPMSEEESSDVEIRNEENNSGLFVTVLEVSSASAAHTG LYTCYYNHTQTEENELEGRHIYIYVDPDPDAFVPLGMTDYLVIVEDDDSAIIPCRITDPETPVTLHNSEGVVPAS YDSRQGFNGTFTVGPYICEATVKGKKFQTIPFNVYALKATSELDLEMEALKTVYKSGETIVVTCAVFNNEVVDL QWTYPGEVKGKGITMLEEIKVPSIKLVYTLTVPPEATVKDSDYECARQATREVKEMKKVTISVHEKGFIEIKPT FSQLEAVNLHEVKHFVVEVRAYPPIRISWLKNNLTIENLITEITDVEKIQEIRYRSKLLKIRAKEEDSGHYTIVAQ NEDAVKSYTFELLTQVPSSILDVDDHHGGSTGGQTVRCTAEGTPLDIEWMICKDIKKCNNETSWTILANNVSN ITEIHSRDRSTVEGRVTFKVEETIAVRCLAKNLLGAENRELKLVAPTLRSEHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only