

Recombinant Human Glucagon,GCG (C-6His)

Catalog No: #AP75302



Package Size: #AP75302-1 10ug #AP75302-2 50ug #AP75302-3 500ug #AP75302-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Glucagon,GCG (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01275
Uniprot	P01275
GeneID	2641;
Calculated MW	18.6kDa
Target Sequence	RSLQDTEEKSRFSASQADPLSDPDQMNEDKRHSQGTFTSDYSKYLDSRRAQDFVQWLMNTKRNRNNAKR HDEFERHAEGTFTSDVSSYLEGQAAKEFIAWLVKGRGRRDFPEEVAIVEELGRRHADGSFSDEMNTILDNLAA RDFINWLIQTKITDRK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,200mM NaCl,1mM DTT,50% Glycerol,pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only