

Recombinant Human ST6 Sialyltransferase 2,ST6GalNAc2 (C-6His)

Catalog No: #AP75304

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75304-1 10ug #AP75304-2 50ug #AP75304-3 500ug #AP75304-4 1mg

Description

| | |
|-----------------|---|
| Product Name | Recombinant Human ST6 Sialyltransferase 2,ST6GalNAc2 (C-6His) |
| Host Species | Human Cells |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | Q9UJ37 |
| Uniprot | Q9UJ37 |
| GeneID | 10610; |
| Calculated MW | 38.84kDa |
| Target Sequence | SAVQRYPGPAAGARDTTSFEAFFQSKASNSWTGKGQACRHLHLAIQRHPHFRGLFNLSIPVLLWGDLFTPAL WDRLSQHKAPYGWRGLSHQVIASLTLNLLNGSESAKLFAPPRDTPPKCIRCAVVGNGGILNGSRQGNIDAHD YVFRNLGAVIKGFERDVGKTSFYGFTVNTMKNSLVSYWNLGFTSVPQGQDLQYIFIPSDIRDYVMLRSAILGV PVPEGLDKGDRPHAYFGPEASASKFKLLHPDFISYLTERRFLKSKLINTHFGDLYMPSTGALMLLTALHTCDQVS AYGFITSNYWKFSDFYFERKMKPLIFYANHDLSEALWRDLHKAGILQLYQR |
| Formulation | Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4. |
| Storage | Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles. |

Note: This product is for in vitro research use only