

## Recombinant Human,Mouse Activin A,INHBA

Catalog No: #AP75315



Package Size: #AP75315-1 10ug #AP75315-2 50ug #AP75315-3 500ug #AP75315-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human,Mouse Activin A,INHBA
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P08476
Uniprot	P08476
GeneID	3624;
Calculated MW	13kDa
Target Sequence	GLECDGKVNICKKQFFVSFKDIGWNDWIIAPSGYHANYCEGECPSHIAGTSGSSLSFHSTVINHYRMRGHSP FANLKSCCVPTKLRPMSMLYYDDGQNIKKDIQNMIVEECGCS
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only