

## Recombinant Mouse Glypican-1, Gpc1 (C-6His)

Catalog No: #AP75320



Package Size: #AP75320-1 10ug #AP75320-2 50ug #AP75320-3 500ug #AP75320-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Glypican-1, Gpc1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9QZF2
Uniprot	Q9QZF2
GeneID	14733;
Calculated MW	56.7kDa
Target Sequence	DPASKSRSCSEVRQIYGAKGFSLSDVPQAEISGEHLRICPQGYTCTSEMEENLANHSRMELESALHDSSRAL QATLATQLHGIDDFQRLLNDSERTLQEAFFGAFGDLYTQNTRAFRDLYAELRLYYRGANLHLEETLAEFWAR LLERLFKQLHPQLLPDDYLDCLGKQAEALRPFGDAPRELRLRATRAFVAARSFVQGLGVASDVVRKVAQVPLA PECSRAIMKLVYCAHCRGVPGARPCPDYCRNVKGLANQADLDAEWRNLLDSMVLITDKFWGPGSGAESVIG GVHVWLAEAINALQDNKDTLTAKVIQACGNPKV/NPHGSGPEEKRRRGKLALQEKPSTGTLEKLVSEAKAQLR DIQDFWISLPGTLCSEKMAMSPASDDRCWNGISKGRYLPEVMGDGLANQINNPEVEVDITKPDMTIRQQIMQL KIMTNRLRGAYGGNDVDFQDASDDGSGSGSGGGCPDDTCGRRVSKKSSSSRTPLTHALPGLSEQEGQKTS VDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only