

Recombinant Mouse C-C Motif Chemokine 24,CCL24,Eotaxin-2

Catalog No: #AP75323

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75323-1 10ug #AP75323-2 50ug #AP75323-3 500ug #AP75323-4 1mg

Description

Product Name	Recombinant Mouse C-C Motif Chemokine 24,CCL24,Eotaxin-2
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9JKC0
Uniprot	Q9JKC0
GeneID	56221;
Calculated MW	10.3kDa
Target Sequence	VTIPSSCCTSFISKKIPENRVVSYQLANGSICPKAGVIFITKKGHKICTDPKLLWVQRHIQKLDKKNQPSKGAKA VRTKFAVQRRRGNSTEV
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only