

Recombinant Mouse Interleukin-10,IL-10

Catalog No: #AP75324



Package Size: #AP75324-1 10ug #AP75324-2 50ug #AP75324-3 500ug #AP75324-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Interleukin-10,IL-10
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P18893
Uniprot	P18893
GeneID	16153;
Calculated MW	18.9kDa
Target Sequence	MSRGQYSREDNNCTHFPVGQSHMLLELRTAFSQVKTFQTKDQLDNILLTDSLMDQDFKGYLGCQALSEMIQF YLVEVMPQAEKHGPEIKEHLNSLGEKLTLMRLRRCHRFLPCENKSKAVEQVKSDFNKLQDQGVYKAMNEF DIFINCIEAYMMIKMKS
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only