

# Recombinant Mouse C-X-C Motif Chemokine 12,CXCL12,SDF-1 $\alpha$

Catalog No: #AP75325

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75325-1 10ug #AP75325-2 50ug #AP75325-3 500ug #AP75325-4 1mg

## Description

Product Name	Recombinant Mouse C-X-C Motif Chemokine 12,CXCL12,SDF-1 $\alpha$
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P40224
Uniprot	P40224
GeneID	20315;
Calculated MW	8.1kDa
Target Sequence	MKPVLSYRCPGRFFESHIANVKHLKILNTPNCALQIVARLKNNNRQVCIDPKLKWIQEYLEKALNK
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only