

Recombinant Human Growth Hormone, GH (Pituitary, 22kD)



Catalog No: #AP75327

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75327-1 10ug #AP75327-2 50ug #AP75327-3 500ug #AP75327-4 1mg

Description

Product Name	Recombinant Human Growth Hormone, GH (Pituitary, 22kD)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01241
Uniprot	P01241
GeneID	2688;
Calculated MW	22.1kDa
Target Sequence	FPTIPLSRLFDNAMLRAHRLHQLAFDITYQEFEEAYIPKEQKYSFLQNPQTSLCFSESIPTPSNREETQQKSNLE LLRISLLLIQSWLEPVQFLRSVFANSLVYGASDSNVYDLLKDLEEGIQTLMGRLEDGSPRTGQIFKQTYSKFDTN SHNDDALLKNYGLLYCFRKMDMKVETFLRIVQCRSVEGSCGF
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only