

# Recombinant Human Dual Specificity Protein Phosphatase 3, DUSP3 (N-6His)



Catalog No: #AP75331

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75331-1 10ug #AP75331-2 50ug #AP75331-3 500ug #AP75331-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Dual Specificity Protein Phosphatase 3, DUSP3 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P51452
Uniprot	P51452
GeneID	1845;
Calculated MW	22.6kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMSGFELSVQDLNDLLSDGSGCYSLPQCNEVTPRIYVGNASVAQDIPKL QKLGITHVLNAAEGRSFMHVNTNANFYKDSGITYLGIKANDTQEFNLSAYFERAADFIDQALAQKNGRVLVHCR EGYSRSPTLVIAYLMMRQKMDVKSALSIVRQNREIGPNDGFLAQLCQLNDRLAKEGKLP
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only