

## Recombinant Mouse Betacellulin,BTC (N-6His)

Catalog No: #AP75333



Package Size: #AP75333-1 10ug #AP75333-2 50ug #AP75333-3 500ug #AP75333-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Betacellulin,BTC (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q543J8
Uniprot	Q543J8
Calculated MW	12.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMDGNTTRTPETNGSLCGAPGENCTGTTTPRQVKVTHFSRCPKQYKHYCIHG RCRFVVD EQTPSCICEKGYFGARCERVDL FYLQQDRGQ
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only