

Recombinant Human Serine Protease HTRA2,HTRA2,Omi (C-6His)



Catalog No: #AP75337

Orders: order@signalwayantibody.com

Package Size: #AP75337-1 10ug #AP75337-2 50ug #AP75337-3 500ug #AP75337-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Serine Protease HTRA2,HTRA2,Omi (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O43464
Uniprot	O43464
GeneID	27429;
Calculated MW	36kDa
Target Sequence	AVPSPPPASPRSQYNFIADVVEKTAPAVVYIEILDRHPFLGREVPISNGSGFVVAADGLIVTNAHVVADRRRVR VRLLSGDTYEAVVTAVDPVADIATLRIQTKEPLPTLPLGRSADVRQGEFVAMGSPFALQNTITSGIVSSAQR ARDLGLPQTNVEYIQTDAIDFGNSGGPLVNL DGEVIGVNTMKVTAGISFAIPSDRLREFLHRGEKKNSSSGIS GSQRRYIGVMMLTLSPSILAEQLREPSFPDVQHGVLIIHKVILGSPAHRAGLRPGDVILAIGE QMVQNAEDVYE AVRTQSQLAVQIRRGRETLTLYVTPEVTEVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only