

Recombinant Human Natural Cytotoxicity Triggering Receptor 1,NCR1,NKp46,CD335 (C-Fc)



Catalog No: #AP75340

Orders: order@signalwayantibody.com

Package Size: #AP75340-1 10ug #AP75340-2 50ug #AP75340-3 500ug #AP75340-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Natural Cytotoxicity Triggering Receptor 1,NCR1,NKp46,CD335 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O76036
Uniprot	O76036
GeneID	9437;
Calculated MW	53.5kDa
Target Sequence	<p>QQQTLPKPFIWAEPHFMPKEKQVTICCGNYGAVEYQLHFEGSLFAVDRPKPPERINKVKFYIPDMNSRMA GQYSCIYRVGELWSEPSNLLDLVTEMYDTPTLSVHPGPEVISGEKVTFYCRLDTATSMFLLLKEGRSSHVQR GYGKVAEFPLGPVTTAHRGTYRCFGSYNNHAWSPFPSEPVKLLVTDIENTSLAPEDPTFPDWTGTYLLTET GLQKDHALWDHTAQNVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVV VDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKT ISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSGDSFFLYS KLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only