

Recombinant Human Parathyroid Hormone 1 Receptor, PTH1R (Gly49, C-6His)

Catalog No: #AP75341

Orders: order@signalwayantibody.com

Package Size: #AP75341-1 10ug #AP75341-2 50ug #AP75341-3 500ug #AP75341-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Parathyroid Hormone 1 Receptor, PTH1R (Gly49, C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q0VGD7
Uniprot	Q0VGD7
Calculated MW	20.2kDa
Target Sequence	YALVDADDVMTKEEQIFLLHRAQAQCGKRLKEVLQRPASIMESDKGWTSASTSGKPRKDKASGKLYPESEED KEAPTGSRYRGRPCLEWDHILCWPLGAPGEVVAVPCPDYIYDFNHKKGHAYRRCDRNGSWELVPGHNRTW ANYSECVKFLTNETREVEVFDRLGMVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only