

Recombinant Mouse Fc γ RII,CD32b,FCGR2 (C-6His)

Catalog No: #AP75342



Package Size: #AP75342-1 10ug #AP75342-2 50ug #AP75342-3 500ug #AP75342-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Fc γ RII,CD32b,FCGR2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P08101
Uniprot	P08101
GeneID	14130;
Calculated MW	21.6kDa
Target Sequence	THDLPKAVVKLEPPWIQVLKEDTVTLTCEGTHNPGNSSTQWFHNGRSIRSQVQASYTFKATVNDSGEYRCQM EQTRLSDPVDLGVISDWLLLQTPQLVFLEGETITLRCHSWRNKLLNRISFFHNEKSVRYHHYSSNFSIPKANHS HSGDYCYCKGSLGRTLHQSKPVTITVQGPKSSRSLPVDHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only