

# Recombinant Human Roundabout Homolog 4,ROBO4 (C-Fc)



Catalog No: #AP75343

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75343-1 10ug #AP75343-2 50ug #AP75343-3 500ug #AP75343-4 1mg

## Description

Product Name	Recombinant Human Roundabout Homolog 4,ROBO4 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8WZ75
Uniprot	Q8WZ75
GeneID	54538;
Calculated MW	74.6kDa
Target Sequence	<p>QDSPPQILVHPQDQLFQGGPARMSCQASGQPPPTIRWLLNGQPLSMVPPDPHLLPDGTLILLQPPARGHA            HDGQALSTDLGVYTCEASNRLGTAVSRGARLSVAVLREDFQIQPRDMVAVVGQFTLECGPPWGHPEPTVS            WWKDGKPLALQPGRHTVSGGSLLMARAESDEGTYMCAVATNSAGHRESRAARVSIQEPQDYTEPVELLA            IQLENVTLNPDPAEGPKPRPAVWLSWKVSGPAAPAQSYTALFRTQTAPGGQGAPWAEELLAGWQSAELGG            LHWGQDYEFKVRPSSGRARGPDSNVLLRLPEKVPSAPPQEVTLKPGNGTVFVSWVPPAENHNGIIRGYQV            WSLGNTSLPPANWTVVGEQTQLEIATHMPGSYCVQAAVTGAGAGEPSRPVCLLLEQAMERATQEPSEHGP            WTLEQLRATLKRVDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVDS            HEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKA            KGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLT            VDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only