

Recombinant Human PRL- 2,PTP4A2 (C-6His)

Catalog No: #AP75348



Package Size: #AP75348-1 10ug #AP75348-2 50ug #AP75348-3 500ug #AP75348-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PRL- 2,PTP4A2 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q12974
Uniprot	Q12974
GeneID	8073;
Calculated MW	20.2kDa
Target Sequence	MNRPAPVEISYENMRFLITHNPTNATLNKFTEELKKYGVTTLVRVCDATYDKAPVEKEGIHVLDWPFDDGAPP PNQIVDDWLNLLKTKFREEPGCCVAVHCVAGLGRAPVLVALALIECGMKYEDAVQFIRQKRRGAFNSKQLLYL EKYRPMRLRFRDRTNGHCCVQLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only