

# Recombinant Human Fibroblast Growth Factor 2,FGF-2,FGFb (Pro143-Ser288)

Catalog No: #AP75350

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75350-1 10ug #AP75350-2 50ug #AP75350-3 500ug #AP75350-4 1mg

## Description

Product Name	Recombinant Human Fibroblast Growth Factor 2,FGF-2,FGFb (Pro143-Ser288)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09038
Uniprot	P09038
GeneID	2247;
Calculated MW	17.3kDa
Target Sequence	PALPEDGGSGAFPPGHFKDPKRLYCKNGGFFLRIHPDGRVDGVREKSDPHIKLQLQAEERGVSIIKGVCANR YLAMKEDGRLLASKCVTDECFERLESNNYNTYRSRKYTSWYVALKRTGQYKLGSKTGPGQKAILFLPMSA KS
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only