

# Recombinant Human Sialic Acid-Binding Ig-Like Lectin 9, Siglec 9, CD329 (C-Fc)



Catalog No: #AP75358

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75358-1 10ug #AP75358-2 50ug #AP75358-3 500ug #AP75358-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Sialic Acid-Binding Ig-Like Lectin 9, Siglec 9, CD329 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y336
Uniprot	Q9Y336
GeneID	27180;
Calculated MW	63.2kDa
Target Sequence	QTSKLLTMQSSVTVQEGLCVHVPCSFSYPSHGWIYPGPVVHGYWFREGANTDQDAPVATNNPARAVWEET RDRFHLLGDPHTENCTLSIRDARRSDAGRYFFRMEKGSIKWNYKHHRLSVNVTALTHRPNILIPGTLESGCPQ NLTCSPWACEQGTTPMISWIGTSVSPLDPSTTRSSVLTLPQPQDHGTSLTCQVTFPGASVTTNKTVHLNVS YPPQNLMTVFQGDGTVSTVLGNGSSLSLPEGQSLRLVCAVDAVDSNPPARLSLSWRGLTLCPSQPSNPGVL ELPWWHLRDEAEFTCRAQNPLGSQQVYLNVSLSKATSGVTQGVDDIEGRMDEPKSCDKTHTCPPCPAPELL GGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLT VLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVE WESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only