

Recombinant Mouse Plexin Domain-Containing Protein 2, PLXDC2 (C-6His)

Catalog No: #AP75360

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75360-1 10ug #AP75360-2 50ug #AP75360-3 500ug #AP75360-4 1mg

Description

Product Name	Recombinant Mouse Plexin Domain-Containing Protein 2, PLXDC2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9DC11
Uniprot	Q9DC11
GeneID	67448;
Calculated MW	48.95kDa
Target Sequence	EPGHHTNDWIYEVTNAFPWNEEGVEVDSQAYNHRWKRNVDPFKAVDTNRASMGQASPESKGFDTLLDDG QDNNTQIEEDTDHNYISRIYGPADSASRDLWVNIDQMEKDKVKIHGILSNTHRQAARVNLSFDFFPYGHFLNE VTVATGGFIYTGVEVHRMLTATQYIAPLMANFDPSVSRNSTVRYFDNGTALVVQWDHVHLQDNYNLGSFTFQ ATLLMDGRIIFGYKEIPVLVTQISSTNHPVKVGLSDFVHVHRIQQIPNVRRTIYEHYRVELQMSKITNISAVEM TPLPTCLQFNGCGPCVSSQIGFNCSWCSKLRQCSGFRHRQDWVDSGCPEEVQSKEKMCKEPTGETSQ TTTTSHTTMQRVLTTRRAVTSQMPTSLPTEDDTKIALHLKDSGASTDDSAEKKGGTLHAVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only