

Recombinant Mouse β -Nerve Growth Factor, β -NGF (Met130-Arg239)

Catalog No: #AP75363

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75363-1 10ug #AP75363-2 50ug #AP75363-3 500ug #AP75363-4 1mg

Description

Product Name	Recombinant Mouse β -Nerve Growth Factor, β -NGF (Met130-Arg239)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P01139
Uniprot	P01139
GeneID	18049;
Calculated MW	12.4kDa
Target Sequence	MGEFSVCDSVSVWVGDKTTATDIKGKEVTVLAEVNINNSVFRQYFFETKCRASNPVESGCRGIDSKHWNSYC TTHTFVKALTTDEKQAAWRFFIRIDTACVCLSRKATR
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM Tris,200mM NaCl,pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only