

Recombinant Human S100 Calcium Binding Protein A8,S100A8,Mrp8 (C-6His)

Catalog No: #AP75364

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75364-1 10ug #AP75364-2 50ug #AP75364-3 500ug #AP75364-4 1mg

Description

Product Name	Recombinant Human S100 Calcium Binding Protein A8,S100A8,Mrp8 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05109
Uniprot	P05109
GeneID	6279;
Calculated MW	11.7kDa
Target Sequence	MLTELEKALNSIIDVYHKYSLIKGNFHAVYRDDLKLLLETECPQYIRKKGADVWFKELDINTDGAVNMQEFLILVI KMGVAAHKKSHEESHKEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM HEPES, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only