

Recombinant Human S100 Calcium Binding Protein A9,S100A9,MRP14



Catalog No: #AP75365

Orders: order@signalwayantibody.com

Package Size: #AP75365-1 10ug #AP75365-2 50ug #AP75365-3 500ug #AP75365-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human S100 Calcium Binding Protein A9,S100A9,MRP14
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P06702
Uniprot	P06702
GeneID	6280;
Calculated MW	13.2kDa
Target Sequence	TCKMSQLERNIETIINTFHQYSVKLGHPDTLNQGEFKELVRKDLQNFLKKENKNEKVIEHIMEDLDTNADKQLSF EEFIMLMARLTWASHEKMHEGDEGPGHHHKPGLGEGTP
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only