

Recombinant Human S100 Calcium Binding Protein A8,A9 Heterodimer (C-6His)



Catalog No: #AP75366

Orders: order@signalwayantibody.com

Package Size: #AP75366-1 10ug #AP75366-2 50ug #AP75366-3 500ug #AP75366-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human S100 Calcium Binding Protein A8,A9 Heterodimer (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05109 & P06702
Uniprot	P05109
GeneID	6279;
Calculated MW	11.7&13.2kDa
Target Sequence	MLTELEKALNSIIDVYHKYSLIKGNFHAVYRDDLKLLLETECPQYIRKKGADVWFKELDINTDGAVNFQEFLILVI KMGVAAHKKSHEESHKEHHHHHH&TCKMSQLERNIETIINTFHQYSVKLGHPDTLNQGEFKELVRKDLQNFLK KENKNEKVIHIMEDLDTNADKQLSFEEFIMLMARLTWASHEKMHEGDEGPGHHHKPGLGEGTP
Formulation	Supplied as a 0.2 µm filtered solution of 20mM HEPES, 10% Glycerol, 150mM NaCl, 2.5mM EDTA, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only