

Recombinant Human, Mouse Fibroblast Growth Factor 8B, FGF-8B



Catalog No: #AP75367

Orders: order@signalwayantibody.com

Package Size: #AP75367-1 10ug #AP75367-2 50ug #AP75367-3 500ug #AP75367-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human, Mouse Fibroblast Growth Factor 8B, FGF-8B
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P55075-3
Uniprot	P55075
GeneID	2253;
Calculated MW	22.5kDa
Target Sequence	MQVTVQSSPNFTQHVREQLVTDQLSRRLIRTYQLYSRTSGKHVQVLANKRINAMAEDGDPFAKLIVETDTFG SRVRVARGAETGLYICMNKKGKLIAKSNGKGDVCFTEIVLENNYTALQNAKYEGWYMAFTRKGRPRKGSKTR QHQREVHFMKRLPRGHHTTEQSLRFEFLNYPFTRSLRGSQRTWAPEPR
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only