

Recombinant Human Galectin-14, LGALS14 (C-6His)

Catalog No: #AP75370



Package Size: #AP75370-1 10ug #AP75370-2 50ug #AP75370-3 500ug #AP75370-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Galectin-14, LGALS14 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8TCE9
Uniprot	Q8TCE9
GeneID	56891;
Calculated MW	17.1kDa
Target Sequence	MSSLPVPYTLPVSLPVGSCVIITGTPILTFVKDPQLEVNIFYTGMDESDIAFQFRLHFGHPAIMNSRVFGIWRYE EKCYLPPFEDGKPFELCIYVRHKEYKVMVNGQRIYNFAHRFPPASVKMLQVLRDISLTRVLISDVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tri, 100mM NaCl, 1mM DTT, 20% glycerol, pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only