

Recombinant Human Ribose-Phosphate Pyrophosphokinase 2, PRPS2 (C-6His)



Catalog No: #AP75371

Orders: order@signalwayantibody.com

Package Size: #AP75371-1 10ug #AP75371-2 50ug #AP75371-3 500ug #AP75371-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ribose-Phosphate Pyrophosphokinase 2, PRPS2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P11908
Uniprot	P11908
GeneID	5634;
Calculated MW	35.8kDa
Target Sequence	PNIVLFGSSHQDLSQRVADRLGLELGKVVTKKFSNQETSVEIGESVRGEDVYIIQSGCGEINDNLMELLIMINA CKIASSSRVTAVIPCFPYARQDKKDKSRAPISAKLVANMLSVAGADHIITMDLHASQIQGFFDIPVDNLYAEPVAVL QWIRENIAEWKNCIIVSPDAGGAKRVTSIADRLNVEFALIHKERKKANEVDRMVLVGDVKDRVAILVDDMADTC GTICHAADKLLSAGATKVYAILTHGIFSGPAISRINNAAFEAVVVTNTIPQEDKMKHCTKIQVIDISMILAEAIRRTH NGESVSYLFSHVPLVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only