

Recombinant Human Galectin 9 (C-6His)

Catalog No: #AP75372



Package Size: #AP75372-1 10ug #AP75372-2 50ug #AP75372-3 500ug #AP75372-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Galectin 9 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q3B8N1-2
Uniprot	Q3B8N1
Calculated MW	36.9kDa
Target Sequence	MAFSGSQAPYLSPAVPFSGTIQGGLQDGLQITVNGTVLSSSGTRFAVNFQTGFSGNDIAFHFNPRFEDGGYVV CNTRQNGSWGPEERKTHMPFQKGMPPDLCFLVQSSDFKVMVNGILFVQYFHRVPPFHRVDTISVNGSVQLSYI SFQPPGVWPANPAPITQTVIHTVQSAPGQMFSTPAIPPMYPHPAYPMPFITLGGLYPSKILLSGTVLPSAQ RFHINLCSGNHIAFHLNPRFDENAVVRNTQIDNSWGSEERSLPRKMPFVRGQSFSVWILCEAHCLKVAVDGQ HLFEYYHRLRNLPINRLEVGGDIQLTHVQTVDDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,5% Trehalose,1mM EDTA,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only