

# Recombinant Human Dipeptidyl-Peptidase 1,DPP1 (C-6His)

Catalog No: #AP75375

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75375-1 10ug #AP75375-2 50ug #AP75375-3 500ug #AP75375-4 1mg

## Description

Product Name	Recombinant Human Dipeptidyl-Peptidase 1,DPP1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16444
Uniprot	P16444
GeneID	1800;
Calculated MW	42.1kDa
Target Sequence	<p>           DFFRDEAERIMRDSVIDGHNDLPWQLLDMFNNRLQDERANLTTLAGTHTNIPKLRAGFVGGQFWSVYTPCD            TQNKDAVRRRTLEQMDVVHRMCRMYPETFLYVTSSAGIRQAFREGKVASLIGVEGGHSIDSSLGVLRALYQLG            MRYLTLTHSCNTPWADNWLVDTGDEPQSQGLSPFGQRVVKELNRLGVLIDLAHVSVATMKATLQLSRAPVI            FSHSSAYSVCASRRNVPDDVLRVLKQTDSLVMVNFYNNYISCTNKANLSQVADHLDHIKEVAGARAVGFGGD            FDGVPRVPEGLEDVSKYPDLIAELLRRNWTEAEVKGALADNLLRVFEAVEQASNLTAPEEEPIPLDQLGGSC            RTHYGYSSVDHHHHHH         </p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only