

# Recombinant Human Airway Trypsin-Like Protease 5,HATL5,TMPRSS11B (C-6His)



Catalog No: #AP75376

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75376-1 10ug #AP75376-2 50ug #AP75376-3 500ug #AP75376-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Airway Trypsin-Like Protease 5,HATL5,TMPRSS11B (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q86T26
Uniprot	Q86T26
GeneID	132724;
Calculated MW	43.1kDa
Target Sequence	LAVEKTYYYQGDFHISGVTYNDNCENAASQASTNLSKDIETKMLNAFQNSSIYKEYVKSEVIKLLPNANGSNVQ LQLKFKFPPAEGVSMRTKIKAKLHQMLKNNMASWNAVPAKIKLMEISKAASEMLTNNCCGRQVANSIITGNKIV NGKSSLEGAWPWQASMQWKGRHYCGASLISSRWLLSAAHCFKAKNNSKDWTVNFGVVVNKPYMTRKVNII FHENYSSPGLHDDIALVQLAEEVSFTEYIRKICLPEAKMKLSENDNVVVTGWGTYMNGSFPVILQEAFKIIDN KICNASYAYSQGFVDSMLCAGFMSGADACQNDSSGGPLAYPDSRNIWHLVGIVSWGDCGCKKNKPGVYTRV TSYRNWITSKTGLVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only