

Recombinant Human Carbonic Anhydrase-Related Protein 11,CA11 (C-6His)



Catalog No: #AP75377

Orders: order@signalwayantibody.com

Package Size: #AP75377-1 10ug #AP75377-2 50ug #AP75377-3 500ug #AP75377-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Carbonic Anhydrase-Related Protein 11,CA11 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O75493
Uniprot	O75493
GeneID	770;
Calculated MW	35.1kDa
Target Sequence	HIGPAPDPEDWWSYKDNLQGNFVPGPPFWGLVNAAWSLCAVGKRQSPVDVELKRVLYDPFLPPLRLSTGGE KLRGTLYNTGRHVSFLPAPRPVNVSGGPLLYSHRLELRLIFGARDGAGSEHQINHQGFSAEVQLIHFNQEL YGNFSAASRGPNGLAILSLFVNVASTSNPFLSRLNDRDITRISYKNDAYFLQDLSLELLFPESFGFITYQGSLST PPCSETVTWILIDRALNITSLQMHSRLRLSQNPSSQIFQSLSGNSRPLQPLAHRALRGNRDPRHPERRCRGPN YRLHVDGVPHGRVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only