

# Recombinant Human Mannoside Acetylglucosaminyltransferase 2, GlcNAc-TII (C-6His)

Catalog No: #AP75378

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75378-1 10ug #AP75378-2 50ug #AP75378-3 500ug #AP75378-4 1mg

## Description

Product Name	Recombinant Human Mannoside Acetylglucosaminyltransferase 2, GlcNAc-TII (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q10469
Uniprot	Q10469
GeneID	4247;
Calculated MW	49.3kDa
Target Sequence	RQRKNEALAPLLDAEPARGAGGRGGDHPVAVGIRRVSNVSAASLVPVAPQPEADNLTLYRSLVYQLNFD QTLRNVDKAGTWAPRELVLVVQVHNRPEYLRLLDLSLRKAQGIDNVLVIFSHDFWSTEINQLIAGVNFPCVLQV FFPFSIQLYPNEFPDPRDCPRDLPKNAALKLGCINAEYPDSFGHYREAKFSQTKHHWWWKLHFVWERVKI LRDYAGLILFLEEDHYLAPDFYHVFKMWWKLQQEQCECDVLSLGTYSASRSFYGMADKVDVKTWKSTEHNM GLALTRNAYQKLIECTDTFCTYDDYNWDWTLQYLTVSCLPKFWKVLVPQIPRIFHAGDCGMHHKTCRPSTQS AQIESLLNNKQYMFPELTISEKFTVVAISPPRKNGGWGDIRDHELCKSYRRLQLDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH8.0 .
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only