

## Recombinant Human GPD1,GDP-C (C-6His)

Catalog No: #AP75379



Package Size: #AP75379-1 10ug #AP75379-2 50ug #AP75379-3 500ug #AP75379-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GPD1,GDP-C (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P21695
Uniprot	P21695
GeneID	2819;
Calculated MW	38.6kDa
Target Sequence	MASKKVCIVGSGNWGSAIAKIVGGNAAQLAQFDPRVTMWWFEEDIGGKKLTEIINTQHENVKYLPGHKLPPNV VAVPDVVQAAEDADILIFVVPHQFIGKICDQLKGHLKANATGISLIKGVDEGPNGLKLISEVIGERLGIPMSVLMG ANIASEVADEKFCETTIGCKDPAQGQLLKELMQTPNFRITVVQEVDTVEICGALKNVVAVGAGFCDGLGFGDN TKAAVIRLGLMEMIAFAKLFCSGPVSSATFLESCGVADLITTCYGGRNRKVAEAFARTGKSIEQLEKELLNGQKL QGPETARELYSILQHKGLVDKFPLFMAVYKVCYEGQPVGFEIHCLQNHPEHMHVDDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only