

Recombinant Human Galectin-3, LGALS3 (C-6His, Human Cells)

Catalog No: #AP75383

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75383-1 10ug #AP75383-2 50ug #AP75383-3 500ug #AP75383-4 1mg

Description

Product Name	Recombinant Human Galectin-3, LGALS3 (C-6His, Human Cells)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P17931
Uniprot	P17931
GeneID	3958;
Calculated MW	27.2kDa
Target Sequence	ADNFSLHDALSGSGNPNPQGWPAGWGNQPAGAGGYPGASYPGAYPGQAPPAYPGQAPPAYPGAPGAYPGAPAGVYPPGPPSGPGAYPSSGQPSATGAYPATGPYGAPAGPLIVPYNLPLPGGVVPRMLITILGTVKPNA NRIALDFQRGNDVAFHFNPRFNENRRVIVCNTKLDNNWGREERQSVFPFESGKPFKIQVLVEPDHFKVAVN DAHLLQYNHRVKKLNEISKLGISGDIDLTSASYTMIVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4, 3mM DTT.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only