

# Recombinant Human GMP Reductase 1,GMPR (C-6His)

Catalog No: #AP75387

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75387-1 10ug #AP75387-2 50ug #AP75387-3 500ug #AP75387-4 1mg

## Description

Product Name	Recombinant Human GMP Reductase 1,GMPR (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P36959
Uniprot	P36959
GeneID	2766;
Calculated MW	38.5kDa
Target Sequence	MPRIDADLKLDKDFKDVLLRPKRSSLKSRAEVDLERTFTFRNSKQTYSGIPIIVANMDTVGTFEMAAVMSQHSMT AIHKHYSLDDWKL FATNHPECLQNAVSSGSGQNDLEKMTSILEAVPQVKFICLDVANGYSEHFVEFVKLVRA KFPEHTIMAGNVVTGEMVEELILSGADIIVGVGPGSVCTTRTKTGVGYPLSAVIECADSAHGLKGHIISDGGC TCPGDVAKAFGAGADFVMLGGMFSGHTECAGEVFERNGRKLKLFYGMSSDTAMNKHAGGVAEYRASEGKT VEVPYKGDVENTILDILGGLRSTCTYVGAALKELSRRTAFIRVTQQHNTVFSVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl,40% glycerol,0.15M NaCl and 1mM DTT, pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only